


# RELIOTRON-based Cathodoluminescence



## **EOS Cathodoluminescence System**


The EOS CL System is a completely integrated optical, digital imaging, fluorescence spectrometer and image analysis/processing package based on a Zeiss Petrographic microscope.

By using the already-proven and widely recognized RELIOTRON chamber/gun/electronic control, this integration allows high resolution CL digital imaging, spectral analysis, and image analysis. Spectral data can be output to Excel, or overlaid from a previous acquisition for direct comparison.

EOS integrators cover the complete system.

## **Light & Electron Optical Chamber Interface**

The vacuum sample holder combines cathodoluminescence observation through a lead glass window with the highest optical efficiency. Microscope objectives can be rotated on the lens turret without moving the chamber. The new electron gun cold cathode source, with a high concentration focus lens coil, provides an ultra stable plasma electron beam focused on the sample.



## **RELIOTRON Specimen Stage**

The CL sample stage is fitted in place of the normal microscope stage. The standard stage has complete field of viewing for a sample up to 70mm x 50mm. X & Y external controls are used to position the specimen.

## **Vacuum System**

The working vacuum for the cold cathode gun is 20 to 100 milli-torr. This pressure is established with a two-stage mechanical rotary pump.



## Specifications

### Stage Standard:

- Slide Drawer for sample load and change
- Concentric X & Y sample position control
- Full view for 70mm x 50mm slides

### Sample Chamber:

- A base window for transmitted light microscopy (TrL)
- Automatic air vent.
- Precision needle gas flow valve
- Hastings pressure gauge
- Chamber cold cathode electron gun assembly
- Industry standard connection for vacuum rotary pump.

### Electronics Console:

- 1-30KV – to 2 mA
- Current stabilized
- Pre-selected gun voltage
- All vacuum interlocked
- Gun overload cut out protection



### EOS Spectrometer

Using a fiber optic interface, and a spectrometer with UV rejecting, high reflectivity, Ag-coated mirrors, provides an enhanced signal from 380-1050 nm.

EOS control software enables spectral capture data output to Excel format.

The graphical overlay history can be stored and used for comparison to live acquisitions.